

AD-A087 995


ENSCO INC SPRINGFIELD VA SIGNAL ANALYSIS SYSTEMS DIV F/G 6/5
ANALOG TO DIGITAL CONVERSION OF ELECTROENCEPHOLOGRAPHIC DATA.(U)
OCT 79 DAAK70-79-M-4337

UNCLASSIFIED

NI

1001
1001

END
DATE
FILMED
9-80
DTIC




MICROCOPY RESOLUTION TEST CHART
NATIONAL BUREAU OF STANDARDS-1963-A


PHOTOGRAPH THIS SHEET

AD A087995

DTIC ACCESSION NUMBER



LEVEL



INVENTORY

ENSCO Inc.,
Springfield VA

Final Rpt., Dtd. 31 Oct. 1979
Contract No. DAAK70-79-M-4337
DOCUMENT IDENTIFICATION

DISTRIBUTION STATEMENT A

Approved for public release
Distribution Unlimited

DISTRIBUTION STATEMENT

ACCESSION FOR	
NTIS	GRA&I <input checked="" type="checkbox"/>
DTIC	TAB <input type="checkbox"/>
UNANNOUNCED	<input type="checkbox"/>
JUSTIFICATION	
Per Ltr. (PL-88, 280-0732)	
BY on file	
DISTRIBUTION	
AVAILABILITY CODES	
DIST	AVAIL AND/OR SPECIAL
A	

DISTRIBUTION STAMP

DTIC
ELECTE
AUG 18 1980
S D D

DATE ACCESSIONED

DATE RECEIVED IN DTIC


PHOTOGRAPH THIS SHEET AND RETURN TO DTIC-DDA-2

AD A087995

ENSEN, INC.

®

80 7 2 094



SIGNAL ANALYSIS SYSTEMS DIVISION

5408A Port Royal Road

Springfield, Virginia 22151

(703) 321-9000

31 October 1979

Reference: Contract DAAK70-79-M-4337

Subject: Final Report on Analog to Digital Conversion of
Electroencephalographic Data

Gentlemen:

This report encompasses the progress made on contract DAAK70-79-M-4337 that commenced 1 August 1979. Services that were to be provided were primarily a level of effort in the analog to digital conversion of 5 Electroencephalographic signal data tapes. We've successfully completed all five tapes, without any major difficulties.

Deliverables were to include the digitized data on 9-track 1/2" magnetic tapes compatible with peripheral equipment located at the MERADCOM computer facility.

After some careful thought the approach taken in digitizing the analog data was to divide the process into four tasks. These tasks are outlined below.

TASK I objectives were to:

1. Rent a recorder capable of reproducing the Electroencephalographic Data. This recorder needed to be outfitted with necessary FM IRIG Low Band, and Direct reproduction electronics, and setup to reproduce the data at a speed of 1 7/8 ips.
2. Purchase a quantity of 10, 1/2" magnetic tapes, allowing a duplicate backup set of digital data tapes to be generated.

TASK II objectives were to:

1. Determine the digitizing specifications to be used for the electroencephalographic data. This step involved the close scrutinization of the analog data with a spectrum analyzer. The analyzer was used in an attempt to determine where most of the

signal power resided in the frequency spectrum over a span of DC to 300Hz. A determination was made that the majority of pertinent information was down below 100Hz, but there was information of possible interest out as far as 200Hz and maybe even 300Hz if sought for vigorously. Therefore initially a small group of events was digitized with the following parameters low pass filter frequency cutoff equal to 300Hz, with a data sample rate of 1024Hz or 1024 points per second, at 1600 BPI.

Upon completion of digitizing this small group of events, it was determined that the next logical step to take would be a magnetic tape dump. This was followed by a times series plot of the data and a power spectrum plot, to verify with a computer, what was observed on the spectrum analyzer. The computer generated times series, and power spectrum plots, truly indicated that the digitization process needed to be altered. The next set was to redo the same events with these new parameters. Low pass filter frequency cutoff, equal to 200Hz, and the sampling rate was dropped to 512Hz, or 512 points per second, at 1600 BPI. These parameters turned out to be ideal for data resolution, and cut in half the amount of digital tape to be used for each analog tape.

2. The next step under this task was to verify that data did exist on the analog tapes. We found it absolutely necessary to redo the log sheets provided with the analog tapes. The decision to redo the log sheets, was made after it was determined, that the time codes supplied for each event with each tape was not at all accurate. This was found to be the case in 90% of the 818 events. The occurrence of the event marker had to be known, and exact to within 1 second. Our objective was to record data in a 10-12 second window, with 5-6 seconds on either side of the event marker. It was found that in some cases time codes logged for event markers was off by as much as 10 seconds and more. It was also discovered that some logged events were never triggered with event markers on the analog tapes. This also proved to be a good reason to check each of the 818 events to be assured of it's identification within the log, and most important, its occurrence.

TASK III objectives were to:

1. Perform mass digitization process, in an effort to complete as many of the 5 analog tapes as possible with funds allocated primarily for this purpose in the contract budget.

Generated deliverables in the form of 5 magnetic digital tapes recorded at 1600BPI on a 9 track drive. Supplied below is information pertaining to digital tape size. Also provided is sufficient documentation to locate and identify specific events on the data tapes. These events will now also be referred to as files.

ANALOG
IDENTIFICATION

NUMBER OF FILES
ON DIGITAL TAPE

Tape 2A	150 Files
Tape 3A	156 Files
Tape 4A	160 Files
Tape 5A	180 Files
Tape 6A	172 Files

Note: First record on file of digital tape 5A is bad. This record has to be skipped to proceed with any processing of electroencepholograhic data.

TAPE # 2A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.	Data
		1	2	3	4	5	START	MAKER	Ses	#	
1	31	TNT	Neutral	EEG-P	EEG-C	--	00:14:00	----	BASELINE		
2	11	TNT	Neutral	EEG-P	EEG-C	--	00:19:53	00:19:58	2-1	A	T
3	11	TNT	Neutral	EEG-P	EEG-C	--	00:20:22	00:20:27	2-1	A	N
4	11	TNT	Neutral	EEG-P	EEG-C	--	00:20:59	00:21:04	2-1	A	T
5	11	TNT	Neutral	EEG-P	EEG-C	--	00:21:25	00:21:30	2-1	A	N
6	11	TNT	Neutral	EEG-P	EEG-C	--	00:21:47	00:21:52	2-1	A	T
7	11	TNT	Neutral	EEG-P	EEG-C	--	00:23:18	00:23:23	2-2	A	T
8	11	TNT	Neutral	EEG-P	EEG-C	--	00:23:43	00:23:48	2-2	A	N
9	11	TNT	Neutral	EEG-P	EEG-C	--	00:24:02	00:24:07	2-2	A	N
10	11	TNT	Neutral	EEG-P	EEG-C	--	00:24:29	00:24:34	2-2	A	T
11	11	TNT	Neutral	EEG-P	EEG-C	--	00:25:05	00:25:10	2-2	A	T
12	11	TNT	Neutral	EEG-P	EEG-C	--	00:26:41	00:26:46	2-3	A	T
13	11	TNT	Neutral	EEG-P	EEG-C	--	00:27:03	00:27:08	2-3	A	N
14	11	TNT	Neutral	EEG-P	EEG-C	--	00:27:42	00:27:47	2-3	A	T
15	11	TNT	Neutral	EEG-P	EEG-C	--	00:28:16	00:28:21	2-3	A	N
16	11	TNT	Neutral	EEG-P	EEG-C	--	00:28:48	00:28:53	2-3	A	T
17	11	TNT	Neutral	EEG-P	EEG-C	--	00:30:20	00:30:25	2-4	A	T
18	11	TNT	Neutral	EEG-P	EEG-C	--	00:30:48	00:30:53	2-4	A	N
19	11	TNT	Neutral	EEG-P	EEG-C	--	00:31:34	00:31:39	2-4	A	T
20	11	TNT	Neutral	EEG-P	EEG-C	--	00:31:59	00:32:04	2-4	A	N
21	11	TNT	Neutral	EEG-P	EEG-C	--	00:31:29	00:33:34	2-5	A	T
22	11	TNT	Neutral	EEG-P	EEG-C	--	00:33:56	00:34:01	2-5	A	N
23	11	TNT	Neutral	EEG-P	EEG-C	--	00:34:21	00:34:26	2-5	A	N
24	11	TNT	Neutral	EEG-P	EEG-C	--	00:34:40	00:34:45	2-5	A	T
25	11	TNT	Neutral	EEGPC	EEG-O-C	--	00:58:07	00:58:12	3-1	A	T
26	11	TNT	Neutral	EEGPC	EEG-O-C	--	00:58:30	00:58:35	3-1	A	N
27	11	TNT	Neutral	EEGPC	EEG-O-C	--	00:58:51	00:58:36	3-1	A	N
28	11	TNT	Neutral	EEGPC	EEG-O-C	--	00:59:12	00:59:17	3-1	A	N
29	11	TNT	Neutral	EEGPC	EEG-O-C	--	00:59:34	00:59:39	3-1	A	T
30	11	TNT	Neutral	EEGPC	EEG-O-C	--	1:01:15	1:01:20	3-2	A	T

TAPE # 2A

DIGITAL INFORMATION						ANALOG INFORMATION					
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec	Rat.	Da
		1	2	3	4	5	START	MAKER	Ses	#	
31	11	TNT	Neutral	EEGP-C	EEGO-C	--	1:01:40	1:01:45	3-2	A	N
32	11	TNT	Neutral	EEGP-C	EEGOC	--	1:01:59	1:02:04	3-2	A	N
33	11	TNT	Neutral	EEGP-C	EEGOC	--	1:02:17	1:02:22	3-2	A	T
34	11	TNT	Neutral	EEGP-C	EEGOC	--	1:02:42	1:02:47	3-2	A	N
35	11	TNT	Neutral	EEGP-C	EEGOC	--	1:03:43	1:03:48	3-3	A	T
36	11	TNT	Neutral	EEGP-C	EEGOC	--	1:04:05	1:04:10	3-3	A	N
37	11	TNT	Neutral	EEGP-C	EEGOC	--	1:04:46	1:04:51	3-3	A	T
38	11	TNT	Neutral	EEGP-C	EEGOC	--	1:05:05	1:05:10	3-3	A	T
39	11	TNT	Neutral	EEGP-C	EEGOC	--	1:05:25	1:05:30	3-3	A	N
40	11	TNT	Neutral	EEGP-C	EEGOC	--	1:07:21	1:07:26	3-4	A	T
41	11	TNT	Neutral	EEGP-C	EEGOC	--	1:07:41	1:07:46	3-4	A	N
42	11	TNT	Neutral	EEGP-C	EEGOC	--	1:07:56	1:08:01	3-4	A	N
43	11	TNT	Neutral	EEGP-C	EEGOC	--	1:08:13	1:08:18	3-4	A	N
44	11	TNT	Neutral	EEGP-C	EEGOC	--	1:08:25	1:08:30	3-4	A	T
45	11	TNT	Neutral	EEGP-C	EEGOC	--	1:09:35	1:09:40	3-5	A	T
46	11	TNT	Neutral	EEGP-C	EEGOC	--	1:09:50	1:09:55	3-5	A	T
47	11	TNT	Neutral	EEGP-C	EEGOC	--	1:10:18	1:10:23	3-5	A	N
48	11	TNT	Neutral	EEGP-C	EEGOC	--	1:10:39	1:10:44	3-5	A	N
49	11	TNT	Neutral	EEGP-C	EEGOC	--	1:10:52	1:10:57	3-5	A	T
50	11	TNT	Neutral	EEGP	EEGO	--	1:33:58	1:34:03	5-1	B	N
51	11	TNT	Neutral	EEGP	EEGO	--	1:34:23	1:34:28	5-1	B	N
52	11	TNT	Neutral	EEGP	EEGO	--	1:34:47	1:34:52	5-1	B	T
53	11	TNT	Neutral	EEGP	EEGO	--	1:35:16	1:35:21	5-1	B	N
54	11	TNT	Neutral	EEGP	EEGO	--	1:35:51	1:35:56	5-1	B	T
55	11	TNT	Neutral	EEGP	EEGO	--	1:36:20	1:36:25	5-1	B	T
56	11	TNT	Neutral	EEGP	EEGO	--	1:39:35	1:39:40	5-2	B	T
57	11	TNT	Neutral	EEGP	EEGO	--	1:40:01	1:40:06	5-2	B	T
58	11	TNT	Neutral	EEGP	EEGO	--	1:40:20	1:40:25	5-2	B	N
59	11	TNT	Neutral	EEGP	EEGO	--	1:40:40	1:40:45	5-2	B	T
60	11	TNT	Neutral	EEGP	EEGO	--	1:42:48	1:42:53	5-3	B	N

TAPE # 2A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec. Rat.	
		1	2	3	4	5	START	MAKER	Ses #	
61	11	TNT	Neutral	EEG P	EEG C	--	1:43:15	1:43:20	5-3	B T
62	11	TNT	Neutral	EEG P	EEG C	--	1:43:37	1:43:42	5-3	B T
63	11	TNT	Neutral	EEG P	EEG C	--	1:44:00	1:44:05	5-3	B N
64	11	TNT	Neutral	EEG P	EEG C	--	1:46:31	1:46:36	5-4	B T
65	11	TNT	Neutral	EEG P	EEG C	--	1:47:28	1:47:33	5-4	B N
66	11	TNT	Neutral	EEG P	EEG C	--	1:47:54	1:47:59	5-4	B T
67	11	TNT	Neutral	EEG P	EEG C	--	1:48:14	1:48:19	5-4	B N
68	11	TNT	Neutral	EEG P	EEG C	--	1:49:31	1:49:36	5-5	B N
69	11	TNT	Neutral	EEG P	EEG C	--	1:50:01	1:50:06	5-5	B T
70	11	TNT	Neutral	EEG P	EEG C	--	1:50:24	1:50:29	5-5	B T
71	11	TNT	Neutral	EEG P	EEG C	--	1:50:53	1:50:58	5-5	B N
72	11	TNT	Neutral	EEG P	EEG C	--	1:51:22	1:51:27	5-5	B T
73	11	TNT	Neutral	EEG P	EEG C	--	1:53:04	1:53:07	5-6	B T
74	11	TNT	Neutral	EEG P	EEG C	--	1:53:27	1:53:32	5-6	B N
75	11	TNT	Neutral	EEG P	EEG C	--	1:53:49	1:53:54	5-6	B N
76	11	TNT	Neutral	EEG P	EEG C	--	1:54:21	1:54:26	5-6	B T
77	11	TNT	Neutral	EEG P	EEG C	--	1:54:44	1:54:49	5-6	B N
78	11	TNT	Neutral	EEG PC	EEGOC	--	2:16:03	2:16:08	6-1	B T
79	11	TNT	Neutral	EEG PC	EEGOC	--	2:16:22	2:16:27	6-1	B T
80	11	TNT	Neutral	EEG PC	EEGOC	--	2:16:44	2:16:49	6-1	B N
81	11	TNT	Neutral	EEG PC	EEGOC	--	2:17:05	2:17:10	6-1	B N
82	11	TNT	Neutral	EEG PC	EEGOC	--	2:17:23	2:17:28	6-1	B T
83	11	TNT	Neutral	EEG PC	EEGOC	--	2:18:23	2:18:28	6-2	B N
84	11	TNT	Neutral	EEG PC	EEGOC	--	2:18:40	2:18:45	6-2	B T
85	11	TNT	Neutral	EEG PC	EEGOC	--	2:18:58	2:19:03	6-2	B T
86	11	TNT	Neutral	EEG PC	EEGOC	--	2:19:19	2:19:24	6-2	B N
87	11	TNT	Neutral	EEG PC	EEGOC	--	2:19:43	2:19:48	6-2	B T
88	11	TNT	Neutral	EEG PC	EEGOC	--	2:21:40	2:21:45	6-3	B N
89	11	TNT	Neutral	EEG PC	EEGOC	--	2:22:08	2:22:13	6-3	B T
90	11	TNT	Neutral	EEG PC	EEGOC	--	2:22:50	2:22:55	6-3	B T

TAPE # 2A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.	Da
		1	2	3	4	5	START	MAKER	Ses	#	
91	11	TNT	Neutral	EEGP-C	EEGOC	---	2:23:11	2:23:16	6-3	B	N
92	11	TNT	Neutral	EEGP-C	EEGOC	---	2:23:30	2:23:35	6-3	B	T
93	11	TNT	Neutral	EEGP-C	EEGOC	---	2:24:42	2:24:47	6-4	B	T
94	11	TNT	Neutral	EEGP-C	EEGOC	---	2:24:59	2:25:04	6-4	B	T
95	11	TNT	Neutral	EEGP-C	EEGOC	---	2:25:21	2:25:26	6-4	B	N
96	11	TNT	Neutral	EEGP-C	EEGOC	---	2:25:50	2:25:55	6-4	B	N
97	11	TNT	Neutral	EEGP-C	EEGOC	---	2:26:06	2:26:11	6-4	B	T
98	11	TNT	Neutral	EEGP-C	EEGOC	---	2:27:17	2:27:22	6-5	B	T
99	11	TNT	Neutral	EEGP-C	EEGOC	---	2:27:36	2:27:41	6-5	B	T
100	11	TNT	Neutral	EEGP-C	EEGOC	---	2:28:03	2:28:08	6-5	B	N
101	11	TNT	Neutral	EEGP-C	EEGOC	---	2:28:25	2:28:30	6-5	B	T
102	11	TNT	Neutral	EEGP-C	EEGOC	---	2:28:45	2:28:50	6-5	B	N
103	11	TNT	Neutral	EEGP-C	EEGOC	---	2:58:21	2:58:26	4-1	A	N
104	11	TNT	Neutral	EEGP-C	EEGOC	---	2:58:37	2:58:42	4-1	A	T
105	11	TNT	Neutral	EEGP-C	EEGOC	---	2:58:53	2:58:58	4-1	A	T
106	11	TNT	Neutral	EEGP-C	EEGOC	---	2:59:20	2:59:25	4-1	A	N
107	11	TNT	Neutral	EEGP-C	EEGOC	---	2:59:32	2:59:37	4-1	A	N
108	11	TNT	Neutral	EEGP-C	EEGOC	---	2:59:48	2:59:53	4-1	A	T
109	11	TNT	Neutral	EEGP-C	EEGOC	---	3:00:47	3:00:52	4-2	A	T
110	11	TNT	Neutral	EEGP-C	EEGOC	---	3:01:04	3:01:09	4-2	A	N
111	11	TNT	Neutral	EEGP-C	EEGOC	---	3:01:30	3:01:35	4-2	A	N
112	11	TNT	Neutral	EEGP-C	EEGOC	---	3:01:52	3:01:57	4-2	A	T
113	11	TNT	Neutral	EEGP-C	EEGOC	---	3:02:09	3:02:14	4-2	A	T
114	11	TNT	Neutral	EEGP-C	EEGOC	---	3:03:37	3:03:42	4-3	A	N
115	11	TNT	Neutral	EEGP-C	EEGOC	---	3:03:57	3:04:02	4-3	A	T
116	11	TNT	Neutral	EEGP-C	EEGOC	---	3:04:18	3:04:23	4-3	A	T
117	11	TNT	Neutral	EEGP-C	EEGOC	---	3:04:36	3:04:41	4-3	A	N
118	11	TNT	Neutral	EEGP-C	EEGOC	---	3:05:00	3:05:05	4-3	A	T
119	11	TNT	Neutral	EEGP-C	EEGOC	---	3:05:16	3:05:21	4-3	A	N
120	11	TNT	Neutral	EEGP-C	EEGOC	---	3:06:20	3:06:25	4-4	A	T

TAPE # 2A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rat.	
		1	2	3	4	5	START	MAKER	Ses #	
121	11	TNT	Neutral	EEG-C	EEG-P	--	3:06:32	3:06:37	4-4	A T
122	11	TNT	Neutral	EEG-C	EEG-P	--	3:06:58	3:07:03	4-4	A N
123	11	TNT	Neutral	EEG-C	EEG-P	--	3:07:14	3:07:19	4-4	A N
124	11	TNT	Neutral	EEG-C	EEG-P	--	3:07:36	3:07:41	4-4	A N
125	11	TNT	Neutral	EEG-C	EEG-P	--	3:07:51	3:07:56	4-4	A T
126	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:29:08	3:29:13	5-1	A T
127	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:29:31	3:29:36	5-1	A N
128	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:29:52	3:29:57	5-1	A N
129	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:30:10	3:30:15	5-1	A T
130	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:30:33	3:30:38	5-1	A T
131	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:32:07	3:32:12	5-2	A T
132	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:32:26	3:12:31	5-2	A T
133	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:32:47	3:32:52	5-2	A N
134	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:33:04	3:33:09	5-2	A T
135	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:33:29	3:33:34	5-2	A T
136	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:52:12	3:35:17	5-3	A N
137	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:36:00	3:36:05	5-3	A T
138	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:36:26	3:36:31	5-3	A T
139	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:36:54	3:36:59	5-3	A N
140	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:37:31	3:37:36	5-3	A N
141	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:39:18	3:39:23	5-4	A T
142	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:39:42	3:39:47	5-4	A T
143	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:40:06	3:40:11	5-4	A N
144	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:40:29	3:40:34	5-4	A T
145	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:41:17	3:41:22	5-4	A N
146	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:42:45	3:42:47	5-5	A N
147	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:42:58	3:43:03	5-5	A N
148	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:43:21	3:43:26	5-5	A T
149	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:43:40	3:43:45	5-5	A T
150	11	TNT	Neutral	EEG P-C	EEGO-C	--	3:44:03	3:44:08	5-5	A N

TAPE # 3A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. #	Stat.
		1	2	3	4	5	START	MAKER	Ses. #	
1	11	TNT	Neutral	EEG OP	EEGP-C	EEG CO	00:24:30	----		Baseline
2	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:29:11	00:29:16	6-1	A T
3	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:29:38	00:29:43	6-1	A T
4	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:30:02	00:30:07	6-1	A N
5	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:30:29	00:30:34	6-1	A N
6	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:30:49	00:30:54	6-1	A T
7	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:31:16	00:31:21	6-1	A N
8	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:33:06	00:33:11	6-2	A N
9	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:33:46	00:33:51	6-2	A T
10	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:34:11	00:34:16	6-2	A T
11	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:34:49	00:34:54	6-2	A N
12	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:35:14	00:35:19	6-2	A T
13	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:35:45	00:35:49	6-2	A N
14	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:37:04	00:37:09	6-3	A T
15	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:37:41	00:37:46	6-3	A T
16	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:38:13	00:38:18	6-3	A N
17	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:38:38	00:38:43	6-3	A T
18	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:39:07	00:39:12	6-3	A N
19	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:39:24	00:39:29	6-3	A N
20	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:41:38	00:41:43	6-4	A N
21	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:42:01	00:42:06	6-4	A T
22	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:42:21	00:42:26	6-4	A T
23	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:42:48	00:42:53	6-4	A N
24	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:43:20	00:43:25	6-4	A T
25	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:43:44	00:43:49	6-4	A N
26	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:44:08	00:44:13	6-4	A N
27	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:46:05	00:46:10	6-5	A T
28	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:46:23	00:46:28	6-5	A T
29	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:46:37	00:46:42	6-5	A N
30	11	TNT	Neutral	EEG OP	EEGP-C	EEGC-O	00:47:08	00:47:13	6-5	A T

TAPE # 3A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.
		1	2	3	4	5	START	MAKER	Ses	#
31	11	TNT	Neutral	EEG OP	EEG PC	EEGC-O	00:47:47	00:47:52	6-5	A N
32	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:07:52	1:07:57	7-1	B T
33	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:08:23	1:08:28	7-1	B N
34	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:08:52	1:08:57	7-1	B N
35	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:09:18	1:09:23	7-1	B T
36	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:09:44	1:09:49	7-1	B T
37	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:10:09	1:10:14	7-1	B N
38	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:11:40	1:11:45	7-2	B T
39	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:12:09	1:12:14	7-2	B N
40	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:12:47	1:12:52	7-2	B N
41	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:13:23	1:13:28	7-2	B T
42	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:13:46	1:13:51	7-2	B T
43	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:15:03	1:15:08	7-3	B T
44	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:15:29	1:15:34	7-3	B N
45	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:16:02	1:16:07	7-3	B T
46	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:16:32	1:16:37	7-3	B T
47	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:16:56	1:17:01	7-3	B N
48	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:17:25	1:17:30	7-3	B T
49	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:18:58	1:19:03	7-4	B N
50	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:19:39	1:19:44	7-4	B T
51	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:20:13	1:20:18	7-4	B N
52	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:21:14	1:21:19	7-4	B T
53	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:22:09	1:22:14	7-4	B N
54	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:22:51	1:22:56	7-4	B T
55	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:24:10	1:24:15	7-5	B N
56	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:24:32	1:24:37	7-5	B N
57	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:24:58	1:25:03	7-5	B T
58	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:25:23	1:25:28	7-5	B T
59	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:26:19	1:26:24	7-5	B T
60	11	TNT	Neutral	EEG P	EEG C	EEG-O	1:26:59	1:27:04	7-5	B N

TAPE # 3A

DIGITAL INFORMATION						ANALOG INFORMATION					
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.	D.
		1	2	3	4	5	START	MAKER	Sec	#	
61	11	TNT	Neutral	EEG C	EEG P	EEG O	1:06:39	1:06:44	1-1	C	N
62	11	TNT	Neutral	EEG C	EEG P	EEG O	1:07:00	1:07:05	1-1	C	N
63	11	TNT	Neutral	EEG C	EEG P	EEG O	1:07:21	1:07:26	1-1	C	T
64	11	TNT	Neutral	EEG C	EEG P	EEG O	1:07:47	1:07:52	1-1	C	N
65	11	TNT	Neutral	EEG C	EEG P	EEG O	1:08:18	1:08:23	1-1	C	T
66	11	TNT	Neutral	EEG C	EEG P	EEG O	1:09:06	1:09:11	1-1	S	N
67	11	TNT	Neutral	EEG C	EEG P	EEG O	1:09:35	1:09:40	1-1	C	T
68	11	TNT	Neutral	EEG C	EEG P	EEG O	1:10:56	1:10:01	1-1	C	T
69	11	TNT	Neutral	EEG C	EEG P	EEG O	1:11:27	1:11:32	1-2	C	T
70	11	TNT	Neutral	EEG C	EEG P	EEG O	1:11:49	1:11:54	1-2	C	T
71	11	TNT	Neutral	EEG C	EEG P	EEG O	1:12:11	1:12:16	1-2	C	N
72	11	TNT	Neutral	EEG C	EEG P	EEG O	1:12:36	1:12:41	1-2	C	N
73	11	TNT	Neutral	EEG C	EEG P	EEG O	1:13:03	1:13:08	1-2	C	N
74	11	TNT	Neutral	EEG C	EEG P	EEG O	1:13:21	1:13:26	1-2	C	T
75	11	TNT	Neutral	EEG C	EEG P	EEG O	1:13:51	1:13:56	1-2	C	T
76	11	TNT	Neutral	EEG C	EEG P	EEG O	1:15:14	1:15:19	1-3	C	T
77	11	TNT	Neutral	EEG C	EEG P	EEG O	1:15:52	1:15:57	1-3	C	N
78	11	TNT	Neutral	EEG C	EEG P	EEG O	1:16:23	1:16:28	1-3	C	N
79	11	TNT	Neutral	EEG C	EEG P	EEG O	1:16:51	1:16:56	1-3	C	T
80	11	TNT	Neutral	EEG C	EEG P	EEG O	1:17:25	1:17:30	1-3	C	N
81	11	TNT	Neutral	EEG C	EEG P	EEG O	1:17:44	1:17:49	1-3	C	N
82	11	TNT	Neutral	EEG C	EEG P	EEG O	1:18:13	1:18:18	1-3	C	T
83	11	TNT	Neutral	EEG C	EEG P	EEG O	1:19:26	1:19:31	1-4	C	T
84	11	TNT	Neutral	EEG C	EEG P	EEG O	1:19:46	1:19:51	1-4	C	T
85	11	TNT	Neutral	EEG C	EEG P	EEG O	1:20:19	1:20:24	1-4	C	T
86	11	TNT	Neutral	EEG C	EEG P	EEG O	1:20:38	1:20:43	1-4	C	N
87	11	TNT	Neutral	EEG C	EEG P	EEG O	1:21:02	1:21:07	1-4	C	N
88	11	TNT	Neutral	EEG C	EEG P	EEG O	1:21:30	1:21:35	1-4	C	T
89	11	TNT	Neutral	EEG C	EEG P	EEG O	1:22:47	1:22:52	1-5	C	N
90	11	TNT	Neutral	EEG C	EEG P	EEG O	1:23:05	1:23:10	1-5	C	N

TAPE # 3A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. #	Stat.	Da
		1	2	3	4	5	START	MAKER			
91	11	TNT	Neutral	EEG C	EEG P	EEG O	1:23:25	1:23:30	1-5	C	T
92	11	TNT	Neutral	EEG C	EEG P	EEG O	1:23:45	1:23:50	1-5	C	T
93	11	TNT	Neutral	EEG C	EEG P	EEG O	1:24:00	1:24:05	1-5	C	N
94	11	TNT	Neutral	EEG C	EEG P	EEG O	1:24:20	1:24:25	1-5	C	T
95	11	TNT	Neutral	EEG C	EEG P	EEG O	1:24:39	1:24:44	1-5	C	T
96	11	TNT	Neutral	EEG C	EEG P	EEG O	1:47:12	1:47:17	2-1	C	T
97	11	TNT	Neutral	EEG C	EEG P	EEG O	1:47:51	1:47:56	2-1	C	N
98	11	TNT	Neutral	EEG C	EEG P	EEG O	1:48:30	1:48:35	2-1	C	N
99	11	TNT	Neutral	EEG C	EEG P	EEG O	1:49:35	1:49:40	2-1	C	T
100	11	TNT	Neutral	EEG C	EEG P	EEG O	1:50:12	1:50:17	2-1	C	T
101	11	TNT	Neutral	EEG C	EEG P	EEG O	1:51:05	1:51:10	2-1	C	N
102	11	TNT	Neutral	EEG P	EEG C	EEG O	1:53:20	1:53:25	2-2	C	N
103	11	TNT	Neutral	EEG P	EEG C	EEG O	1:53:56	1:54:01	2-2	C	N
104	11	TNT	Neutral	EEG P	EEG C	EEG O	1:54:30	1:54:35	2-2	C	T
105	11	TNT	Neutral	EEG P	EEG C	EEG O	1:55:18	1:55:23	2-2	C	T
106	11	TNT	Neutral	EEG P	EEG C	EEG O	1:56:45	1:56:50	2-2	C	N
107	11	TNT	Neutral	EEG P	EEG C	EEG O	1:57:15	1:57:20	2-2	C	T
108	11	TNT	Neutral	EEG P	EEG C	EEG O	1:59:25	1:59:30	2-2	C	T
109	11	TNT	Neutral	EEG P	EEG C	EEG O	1:59:51	1:59:56	2-3	C	N
110	11	TNT	Neutral	EEG P	EEG C	EEG O	2:00:13	2:00:18	2-3	C	N
111	11	TNT	Neutral	EEG P	EEG C	EEG O	2:00:36	2:00:41	2-3	C	T
112	11	TNT	Neutral	EEG P	EEG C	EEG O	2:00:57	2:01:02	2-3	C	N
113	11	TNT	Neutral	EEG P	EEG C	EEG O	2:01:27	2:01:32	2-3	C	T
114	11	TNT	Neutral	EEG P	EEG C	EEG O	2:03:09	2:03:14	2-4	C	T
115	11	TNT	Neutral	EEG P	EEG C	EEG O	2:03:34	2:03:39	2-4	C	N
116	11	TNT	Neutral	EEG P	EEG C	EEG O	2:04:00	2:04:05	2-4	C	N
117	11	TNT	Neutral	EEG P	EEG C	EEG O	2:04:23	2:04:28	2-4	C	T
118	11	TNT	Neutral	EEG P	EEG C	EEG O	2:04:56	2:05:01	2-4	C	T
119	11	TNT	Neutral	EEG P	EEG C	EEG O	2:05:34	2:05:39	2-4	C	N
120	11	TNT	Neutral	EEG P	EEG C	EEG O	2:07:34	2:07:39	2-5	C	N

TAPE # 3A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.
		1	2	3	4	5	START	MAKER	Ses	#
121	11	TNT	Neutral	EEG P	EEG C	EEG O	2:08:07	2:08:12	2-5	C T
122	11	TNT	Neutral	EEG P	EEG C	EEG O	2:08:58	2:09:03	2-5	C N
123	11	TNT	Neutral	EEG P	EEG C	EEG O	2:09:45	2:09:50	2-5	C T
124	11	TNT	Neutral	EEG P	EEG C	EEG O	2:10:22	2:10:27	2-5	C N
125	11	TNT	Neutral	EEG P	EEG C	EEG O	2:10:51	2:10:56	2-5	C T
126	11	TNT	Neutral	EEG P	EEG C	EEG O	2:11:25	2:11:30	2-5	C T
127	11	TNT	Neutral	EEG P	EEG O	EEG C	2:43:01	2:43:06	3-1	C T
128	11	TNT	Neutral	EEG P	EEG O	EEG C	2:43:47	2:43:52	3-1	C N
129	11	TNT	Neutral	EEG P	EEG O	EEG C	2:44:27	2:44:32	3-1	C T
130	11	TNT	Neutral	EEG P	EEG O	EEG C	2:45:02	2:45:07	3-1	C T
131	11	TNT	Neutral	EEG P	EEG O	EEG C	2:45:37	2:45:42	3-1	C N
132	11	TNT	Neutral	EEG P	EEG O	EEG C	2:46:12	2:46:17	3-1	C T
133	11	TNT	Neutral	EEG P	EEG O	EEG C	2:47:16	2:47:21	3-2	C N
134	11	TNT	Neutral	EEG P	EEG O	EEG C	2:47:41	2:47:46	3-2	C N
135	11	TNT	Neutral	EEG P	EEG O	EEG C	2:48:13	2:48:18	3-2	C T
136	11	TNT	Neutral	EEG P	EEG O	EEG C	2:48:46	2:48:51	3-2	C T
137	11	TNT	Neutral	EEG P	EEG O	EEG C	2:49:18	2:49:23	3-2	C T
138	11	TNT	Neutral	EEG P	EEG O	EEG C	2:49:37	2:49:42	3-2	C N
139	11	TNT	Neutral	EEG P	EEG O	EEG C	2:50:50	2:50:55	3-3	C T
140	11	TNT	Neutral	EEG P	EEG O	EEG C	2:51:25	2:51:30	3-3	C T
141	11	TNT	Neutral	EEG P	EEG O	EEG C	2:52:02	2:52:07	3-3	C N
142	11	TNT	Neutral	EEG P	EEG O	EEG C	2:52:32	2:52:37	3-3	C N
143	11	TNT	Neutral	EEG P	EEG O	EEG C	2:52:59	2:53:04	3-3	C N
144	11	TNT	Neutral	EEG P	EEG O	EEG C	2:53:28	2:53:33	3-3	C T
145	11	TNT	Neutral	EEG P	EEG O	EEG C	2:54:32	2:54:37	3-4	C T
146	11	TNT	Neutral	EEG P	EEG O	EEG C	2:55:00	2:55:05	3-4	C T
147	11	TNT	Neutral	EEG P	EEG O	EEG C	2:55:25	2:55:30	3-4	C N
148	11	TNT	Neutral	EEG P	EEG O	EEG C	2:55:46	2:55:51	3-4	C N
149	11	TNT	Neutral	EEG P	EEG O	EEG C	2:56:00	2:56:05	3-4	C T
150	11	TNT	Neutral	EEG P	EEG O	EEG C	2:56:18	2:56:23	3-4	C N

TAPE # 3A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Stat.
		1	2	3	4	5	START	MAKER	Ses #	
151	11	TNT	Neutral	EEG P	EEG O	EEG C	2:57:27	2:57:32	3-5 C	N
152	11	TNT	Neutral	EEG P	EEG O	EEG C	2:57:57	2:58:02	3-5 C	T
153	11	TNT	Neutral	EEG P	EEG O	EEG C	2:58:13	2:58:18	3-5 C	T
154	11	TNT	Neutral	EEG P	EEG O	EEG C	2:58:29	2:58:34	3-5 C	N
155	11	TNT	Neutral	EEG P	EEG O	EEG C	2:58:48	2:58:53	3-5 C	T
156	11	TNT	Neutral	EEG P	EEG O	EEG C	2:59:24	2:59:29	3-5 C	N

TAPE # 4A

DIGITAL INFORMATION						ANALOG INFORMATION			
File # of	Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec'd at. Da
		1	2	3	4	5	START	MAKER	Ses #
1	11	TNT	Neutral	EEG O	EEG P	EEG C	00:24:42	-----	Baseline
2	11	TNT	Neutral	EEG O	EEG P	EEG C	00:32:19	00:32:24	1-P 98 T
3	11	TNT	Neutral	EEG O	EEG P	EEG C	00:37:50	00:32:55	1-P 98 N
4	11	TNT	Neutral	EEG O	EEG P	EEG C	00:33:17	00:33:22	1-P 98 N
5	11	TNT	Neutral	EEG O	EEG P	EEG C	00:33:44	00:33:49	1-P 98 T
6	11	TNT	Neutral	EEG O	EEG P	EEG C	00:34:23	00:34:28	1-P 98 N
7	11	TNT	Neutral	EEG O	EEG P	EEG C	00:35:17	00:35:22	1-P 98 T
8	11	TNT	Neutral	EEG O	EEG P	EEG C	00:36:00	00:36:05	1-P 98 N
9	11	TNT	Neutral	EEG O	EEG P	EEG C	00:36:28	00:36:33	1-P 98 T
10	11	TNT	Neutral	EEG O	EEG P	EEG C	00:37:05	00:37:10	1-P 98 N
11	11	TNT	Neutral	EEG O	EEG P	EEG C	00:37:37	00:37:42	1-P 98 T
12	11	TNT	Neutral	EEG O	EEG P	EEG C	00:38:03	00:38:08	1-P 98 N
13	11	TNT	Neutral	EEG O	EEG P	EEG C	00:38:35	00:38:40	1-P 98 T
14	11	TNT	Neutral	EEG O	EEG P	EEG C	00:39:11	00:39:16	1-P 98 T
15	11	TNT	Neutral	EEG O	EEG P	EEG C	00:39:37	00:39:42	1-P 98 N
16	11	TNT	Neutral	EEG O	EEG P	EEG C	00:40:22	00:40:27	1-P 98 T
17	11	TNT	Neutral	EEG O	EEG P	EEG C	00:41:01	00:41:06	1-P 98 T
18	11	TNT	Neutral	EEG O	EEG P	EEG C	00:41:50	00:41:55	1-P 98 N
19	11	TNT	Neutral	EEG O	EEG P	EEG C	00:42:26	00:42:31	1-P 98 N
20	11	TNT	Neutral	EEG O	EEG P	EEG C	00:42:55	00:43:00	1-P 98 T
21	11	TNT	Neutral	EEG O	EEG P	EEG C	00:43:15	00:43:20	1-P 98 T
22	11	TNT	Neutral	EEG O	EEG P	EEG C	00:43:32	00:43:37	1-P 98 N
23	11	TNT	Neutral	EEG O	EEG P	EEG C	00:43:57	00:44:02	1-P 98 T
24	11	TNT	Neutral	EEG O	EEG P	EEG C	00:44:31	00:44:36	1-P 98 N
25	11	TNT	Neutral	EEG O	EEG P	EEG C	00:45:05	00:45:10	1-P 98 N
26	11	TNT	Neutral	EEG O	EEG P	EEG C	00:45:39	00:45:44	1-P 98 T
27	11	TNT	Neutral	EEG O	EEG P	EEG C	00:46:01	00:46:06	1-P 98 N
28	11	TNT	Neutral	EEG O	EEG P	EEG C	00:46:27	00:46:32	1-P 98 N
29	11	TNT	Neutral	EEG O	EEG P	EEG C	00:46:46	00:46:51	1-P 98 T
30	11	TNT	Neutral	EEG O	EEG P	EEG C	00:50:11	00:50:16	1-1 98 N

TAPE # 4A

DIGITAL INFORMATION							ANALOG INFORMATION		
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.
		1	2	3	4	5	START	MARKER	Sess #
31	11	TNT	Neutra	EEG O	EEG P	EEG C	00:50:47	00:50:52	1-1 98 T
32	11	TNT	Neutra	EEG O	EEG P	EEG C	00:51:23	00:51:28	1-1 98 N
33	11	TNT	Neutra	EEG O	EEG P	EEG C	00:52:04	00:52:09	1-1 98 T
34	11	TNT	Neutra	EEG O	EEG P	EEG C	00:52:41	00:52:46	1-1 98 T
35	11	TNT	Neutra	EEG O	EEG P	EEG C	00:53:07	00:53:12	1-1 98 N
36	11	TNT	Neutra	EEG O	EEG P	EEG C	00:55:27	00:55:32	1-2 98 N
37	11	TNT	Neutra	EEG O	EEG P	EEG C	00:56:18	00:56:23	1-2 98 T
38	11	TNT	Neutra	EEG O	EEG P	EEG C	00:56:49	00:56:54	1-2 98 N
39	11	TNT	Neutra	EEG O	EEG P	EEG C	00:57:32	00:57:37	1-2 98 T
40	11	TNT	Neutra	EEG O	EEG P	EEG C	00:58:11	00:58:16	1-2 98 N
41	11	TNT	Neutra	EEG O	EEG P	EEG C	00:58:43	00:58:48	1-2 98 T
42	11	TNT	Neutra	EEG O	EEG P	EEG C	1:01:12	1:01:17	1-3 98 N
43	11	TNT	Neutra	EEG O	EEG P	EEG C	1:01:53	1:01:58	1-3 98 T
44	11	TNT	Neutra	EEG O	EEG P	EEG C	1:02:27	1:02:32	1-3 98 N
45	11	TNT	Neutra	EEG O	EEG P	EEG C	1:02:46	1:02:51	1-3 98 N
46	11	TNT	Neutra	EEG O	EEG P	EEG C	1:03:03	1:03:08	1-3 98 T
47	11	TNT	Neutra	EEG O	EEG P	EEG C	1:03:20	1:03:25	1-3 98 N
48	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:27:29	1:27:34	2-1 98 T
49	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:28:06	1:28:11	2-1 98 N
50	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:28:41	1:28:46	2-1 98 T
51	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:29:12	1:29:17	2-1 98 N
52	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:29:38	1:29:43	2-1 98 T
53	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:30:04	1:30:09	2-1 98 T
54	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:31:58	1:32:03	2-2 98 N
55	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:32:16	1:32:21	2-2 98 N
56	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:32:59	1:33:04	2-2 98 T
57	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:33:37	1:33:42	2-2 98 T
58	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:33:57	1:34:02	2-2 98 N
59	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:34:20	1:34:25	2-2 98 T
60	11	TNT	Neutra	EEG P-C	EEG C-O	EEGP-O	1:37:46	1:37:51	2-3 98 T

TAPE 1 4A

DIGITAL INFORMATION						ANALOG INFORMATION					
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec. Ses.	Rat. #	Data
		1	2	3	4	5	START	MAKER			
61	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:38:21	1:38:26	2-3	98	N
62	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:38:44	1:38:49	2-3	98	N
63	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:39:04	1:39:09	2-3	98	T
64	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:39:26	1:39:31	2-3	98	T
65	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:39:47	1:39:52	2-3	98	N
66	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:41:38	1:41:43	2-4	98	T
67	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:42:02	1:42:07	2-4	98	N
68	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:42:22	1:42:27	2-4	98	T
69	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:42:41	1:42:46	2-4	98	N
70	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:42:58	1:43:03	2-4	98	N
71	11	TNT	Neutral	EEG P-C	EEG CO	EEGP-O	1:43:20	1:43:25	2-4	98	T
72	11	TNT	Neutral	EEG C	EEG O	EEG P	2:14:59	2:15:04	4-1	C	N
73	11	TNT	Neutral	EEG C	EEG O	EEG P	2:16:36	2:16:41	4-1	C	N
74	11	TNT	Neutral	EEG C	EEG O	EEG P	2:16:57	2:17:02	4-1	C	T
75	11	TNT	Neutral	EEG C	EEG O	EEG P	2:17:19	2:17:24	4-1	C	T
76	11	TNT	Neutral	EEG C	EEG O	EEG P	2:17:53	2:17:58	4-1	C	N
77	11	TNT	Neutral	EEG C	EEG O	EEG P	2:18:52	2:18:16	4-1	C	T
78	11	TNT	Neutral	EEG C	EEG O	EEG P	2:19:32	2:19:37	4-2	C	T
79	11	TNT	Neutral	EEG C	EEG O	EEG P	2:20:13	2:20:18	4-2	C	N
80	11	TNT	Neutral	EEG C	EEG O	EEG P	2:20:57	2:21:02	4-2	C	T
81	11	TNT	Neutral	EEG C	EEG O	EEG P	2:21:28	2:21:33	4-2	C	N
82	11	TNT	Neutral	EEG C	EEG O	EEG P	2:22:08	2:22:13	4-2	C	T
83	11	TNT	Neutral	EEG C	EEG O	EEG P	2:22:52	2:22:57	4-2	C	N
84	11	TNT	Neutral	EEG C	EEG O	EEG P	2:25:00	2:25:05	4-3	C	N
85	11	TNT	Neutral	EEG C	EEG O	EEG P	2:25:42	2:25:47	4-3	C	T
86	11	TNT	Neutral	EEG C	EEG O	EEG P	2:26:10	2:26:15	4-3	C	T
87	11	TNT	Neutral	EEG C	EEG O	EEG P	2:26:32	2:26:37	4-3	C	N
88	11	TNT	Neutral	EEG C	EEG O	EEG P	2:26:59	2:27:04	4-3	C	T
89	11	TNT	Neutral	EEG C	EEG O	EEG P	2:27:18	2:27:23	4-3	C	T
90	11	TNT	Neutral	EEG C	EEG O	EEG P	2:28:51	2:28:56	4-4	C	N

TAPE # 4A

DIGITAL INFORMATION						ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.
		1	2	3	4	5	START	MAKER	Ses. #
91	11	TNT	Neutral	EEG C	EEG O	EEG P	2:29:25	2:29:30	4-4 C N
92	11	TNT	Neutral	EEG C	EEG O	EEG P	2:29:47	2:29:52	4-4 C T
93	11	TNT	Neutral	EEG C	EEG O	EEG P	2:30:16	2:30:21	4-4 C N
94	11	TNT	Neutral	EEG C	EEG O	EEG P	2:30:59	2:31:04	4-4 C T
95	11	TNT	Neutral	EEG C	EEG O	EEG P	2:31:32	2:31:37	4-4 C N
96	11	TNT	Neutral	EEG C	EEG O	EEG P	2:33:17	2:33:22	4-5 C T
97	11	TNT	Neutral	EEG C	EEG O	EEG P	2:33:57	2:34:02	4-5 C N
98	11	TNT	Neutral	EEG C	EEG O	EEG P	2:34:17	2:34:22	4-5 C N
99	11	TNT	Neutral	EEG C	EEG O	EEG P	2:34:49	2:34:54	4-5 C T
100	11	TNT	Neutral	EEG C	EEG O	EEG P	2:35:29	2:35:34	4-5 C N
101	11	TNT	Neutral	EEG C	EEG O	EEG P	2:36:29	2:36:34	4-5 C T
102	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:11:33	3:11:38	3-1 98 T
103	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:12:19	3:12:24	3-1 98 N
104	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:12:54	3:12:59	3-1 98 N
105	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:13:18	3:13:23	3-1 98 T
106	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:13:47	3:13:52	3-1 98 T
107	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:14:21	3:14:26	3-1 98 T
108	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:15:46	3:15:51	3-2 98 T
109	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:16:12	3:16:17	3-2 98 N
110	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:16:48	3:16:53	3-2 98 T
111	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:17:22	3:17:27	3-2 98 N
112	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:18:10	3:18:15	3-2 98 N
113	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:18:43	3:18:48	3-2 98 T
114	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:19:58	3:20:03	3-3 98 N
115	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:20:23	3:20:28	3-3 98 N
116	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:20:44	3:20:49	3-3 98 T
117	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:21:02	3:21:07	3-3 98 T
118	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:21:27	3:21:32	3-3 98 N
119	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:21:54	3:21:59	3-3 98 T
120	11	TNT	Neutral	EEG CO	EEGPO	EEG PC	3:23:02	3:23:07	3-4 98 N

TAPE 4A

DIGITAL INFORMATION						ANALOG INFORMATION					
File #	of	CONTENTS IN-CHANNEL					TIME CODE		Rec	Rat.	D.
		1	2	3	4	5	START	MAKER	Ses	#	
121	11	TNT	Neutral	EEG CO	EEG PO	EEGPC	3:23:22	3:23:27	3-4	98	N
122	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:23:54	3:23:59	3-4	98	T
123	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:24:16	3:24:21	3-4	98	T
124	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:24:39	3:24:44	3-4	98	N
125	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:25:02	3:25:07	3-4	98	T
126	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:26:29	3:26:34	3-5	98	T
127	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:27:00	3:27:05	3-5	98	T
128	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:27:20	3:27:25	3-5	98	N
129	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:27:39	3:27:44	3-5	98	N
130	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:27:57	3:28:02	3-5	98	T
131	11	TNT	Neutral	EEG CO	EEG PO	EEG PC	3:28:18	3:28:23	3-5	98	N
132	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:50:56	3:51:01	5-1	C	T
133	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:51:36	3:51:41	5-1	C	T
134	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:52:13	3:52:18	5-1	C	N
135	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:52:41	3:52:46	5-1	C	N
136	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:53:26	3:53:31	5-1	C	T
137	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:53:52	3:53:57	5-1	C	N
138	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:55:24	3:55:29	5-2	C	N
139	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:55:52	3:55:57	5-2	C	T
140	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:56:17	3:56:23	5-2	C	T
141	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:56:51	3:56:56	5-2	C	N
142	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:57:25	3:57:30	5-2	C	T
143	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:57:53	3:57:58	5-2	C	N
144	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:59:13	3:59:18	5-2	C	T
145	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	3:59:42	3:59:47	5-3	C	N
146	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:00:18	4:00:23	5-3	C	T
147	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:00:42	4:00:47	5-3	C	N
148	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:01:29	4:01:34	5-3	C	T
149	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:02:57	4:03:02	5-4	C	N
150	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:03:28	4:03:33	5-4	C	T

TAPE # 4A

DIGITAL INFORMATION						ANALOG INFORMATION					
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec	Rat.	
		1	2	3	4	5	START	MAKER	Ses	#	
151	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:04:10	4:04:15	5-4	C	T
152	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:04:39	4:04:44	5-4	C	N
153	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:05:04	4:05:09	5-4	C	T
154	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:05:28	4:05:33	5-4	C	N
155	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:06:58	4:07:03	5-5	C	T
156	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:07:27	4:07:32	5-5	C	N
157	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:07:49	4:07:54	5-5	C	T
158	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:08:12	4:08:17	5-5	C	N
159	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:08:36	4:08:41	5-5	C	T
160	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	4:09:01	4:09:06	5-5	C	N

TAPE # 5A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec. #	Rat.	Dir.
		1	2	3	4	5	START	MAKER			
1	11	TNT	Neutral	EEG O	EEG P	EEG C	00:27:00	N/A	1	G	B
2	11	TNT	Neutral	EEG O	EEG P	EEG C	00:30:50	00:30:55	1-1	G	T
3	11	TNT	Neutral	EEG O	EEG P	EEG C	00:31:42	00:31:47	1-1	G	N
4	11	TNT	Neutral	EEG O	EEG P	EEG C	00:32:25	00:32:30	1-1	G	N
5	11	TNT	Neutral	EEG O	EEG P	EEG C	00:33:10	00:33:15	1-1	G	T
6	11	TNT	Neutral	EEG O	EEG P	EEG C	00:33:55	00:34:00	1-1	G	N
7	11	TNT	Neutral	EEG O	EEG P	EEG C	00:34:40	00:34:45	1-1	G	T
8	11	TNT	Neutral	EEG O	EEG P	EEG C	00:35:19	00:35:24	1-2	G	T
9	11	TNT	Neutral	EEG O	EEG P	EEG C	00:35:55	00:36:00	1-2	G	N
10	11	TNT	Neutral	EEG O	EEG P	EEG C	00:36:58	00:37:03	1-2	G	T
11	11	TNT	Neutral	EEG O	EEG P	EEG C	00:37:47	00:37:52	1-2	G	N
12	11	TNT	Neutral	EEG O	EEG P	EEG C	00:38:35	00:38:40	1-2	G	T
13	11	TNT	Neutral	EEG O	EEG P	EEG C	00:39:16	00:39:21	1-2	G	N
14	11	TNT	Neutral	EEG O	EEG P	EEG C	-----	-----	--	--	-
15	11	TNT	Neutral	EEG O	EEG P	EEG C	00:40:05	00:40:10	1-3	G	N
16	11	TNT	Neutral	EEG O	EEG P	EEG C	00:40:34	00:40:39	1-3	G	N
17	11	TNT	Neutral	EEG O	EEG P	EEG C	00:40:53	00:40:58	1-3	G	T
18	11	TNT	Neutral	EEG O	EEG P	EEG C	00:41:10	00:41:15	1-3	G	T
19	11	TNT	Neutral	EEG O	EEG P	EEG C	00:41:29	00:41:34	1-3	G	N
20	11	TNT	Neutral	EEG O	EEG P	EEG C	00:41:54	00:41:59	1-3	G	T
21	11	TNT	Neutral	EEG O	EEG P	EEG C	00:42:15	00:42:20	1-4	G	T
22	11	TNT	Neutral	EEG O	EEG P	EEG C	00:42:30	00:42:35	1-4	G	T
23	11	TNT	Neutral	EEG O	EEG P	EEG C	00:43:00	00:43:05	1-4	G	N
24	11	TNT	Neutral	EEG O	EEG P	EEG C	00:43:26	00:43:31	1-4	G	N
25	11	TNT	Neutral	EEG O	EEG P	EEG C	00:44:00	00:44:05	1-4	G	N
26	11	TNT	Neutral	EEG O	EEG P	EEG C	00:44:55	00:45:00	1-4	G	T
27	11	TNT	Neutral	EEG O	EEG P	EEG C	00:45:40	00:45:45	1-4	G	T
28	11	TNT	Neutral	EEG O	EEG P	EEG C	00:46:00	00:46:05	1-4	G	T
29	11	TNT	Neutral	EEG O	EEG P	EEG C	00:46:21	00:46:26	1-4	G	N
30	11	TNT	Neutral	EEG O	EEG P	EEG C	00:47:09	00:47:14	1-4	G	N

TAPE # 5A

DIGITAL INFORMATION						ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec. Rat.
		1	2	3	4	5	START	MAKER	Ses. #
31	11	TNT	Neutral	EEG O	EEG P	EEG C	00:47:42	00:47:47	1-5P G T
32	11	TNT	Neutral	EEG O	EEG P	EEG C	00:48:07	00:48:12	1-5P G N
33	18	TNT	Neutral	EEG O	EEG P	EEG C	00:50:00	00:50:05	1-1 G T
34	11	TNT	Neutral	EEG O	EEG P	EEG C	00:50:28	00:50:33	1-1 G T
35	11	TNT	Neutral	EEG O	EEG P	EEG C	00:51:05	00:51:10	1-1 G N
36	11	TNT	Neutral	EEG O	EEG P	EEG C	00:52:00	00:52:05	1-1 G N
37	11	TNT	Neutral	EEG O	EEG P	EEG C	00:52:37	00:52:42	1-1 G T
38	11	TNT	Neutral	EEG O	EEG P	EEG C	00:53:09	00:53:14	1-1 G N
39	11	TNT	Neutral	EEG O	EEG P	EEG C	00:54:35	00:54:40	1-2 G N
40	11	TNT	Neutral	EEG O	EEG P	EEG C	00:54:55	00:55:00	1-2 G N
41	11	TNT	Neutral	EEG O	EEG P	EEG C	00:55:16	00:55:21	1-2 G T
42	11	TNT	Neutral	EEG O	EEG P	EEG C	00:55:39	00:55:44	1-2 G T
43	11	TNT	Neutral	EEG O	EEG P	EEG C	00:56:16	00:56:21	1-2 G N
44	11	TNT	Neutral	EEG O	EEG P	EEG C	00:56:38	00:56:43	1-2 G T
45	11	TNT	Neutral	EEG O	EEG P	EEG C	00:57:55	00:58:00	1-3 G T
46	11	TNT	Neutral	EEG O	EEG P	EEG C	00:58:20	00:58:25	1-3 G T
47	11	TNT	Neutral	EEG O	EEG P	EEG C	00:59:00	00:59:05	1-3 G N
48	11	TNT	Neutral	EEG O	EEG P	EEG C	00:59:36	00:59:41	1-3 G N
49	11	TNT	Neutral	EEG O	EEG P	EEG C	1:00:08	1:00:13	1-3 G T
50	11	TNT	Neutral	EEG O	EEG P	EEG C	1:00:44	1:00:49	1-3 G N
51	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:32:17	1:32:22	2-1 G T
52	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:32:51	1:32:56	2-1 G N
53	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:33:29	1:33:34	2-1 G T
54	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:34:23	1:34:28	2-1 G N
55	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:35:06	1:35:11	2-1 G T
56	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:35:42	1:35:47	2-1 G N
57	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:37:16	1:37:21	2-1 G N
58	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:37:45	1:37:50	2-1 G T
59	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:38:37	1:38:42	2-1 G N
60	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:39:58	1:40:03	2-1 G T

TAPE #5A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Reclat.		
		1	2	3	4	5	START	MAKER	Ses	#	
61	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:40:22	1:40:27	2-2	G	N
62	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:40:45	1:40:50	2-2	G	T
63	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:43:18	1:43:23	2-3	G	T
64	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:43:48	1:43:53	2-3	G	T
65	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:44:19	1:44:24	2-3	G	N
66	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:44:42	1:44:47	2-3	G	N
67	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:45:10	1:45:15	2-3	G	T
68	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:45:36	1:45:41	2-3	G	N
69	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:47:07	1:47:12	2-4	G	T
70	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:47:41	1:47:46	2-4	G	N
71	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:48:16	1:48:21	2-4	G	T
72	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:49:00	1:49:05	2-4	G	N
73	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:49:32	1:49:37	2-4	G	T
74	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:49:54	1:49:59	2-4	G	N
75	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:51:26	1:51:31	2-5	G	T
76	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:52:09	1:52:14	2-5	G	N
77	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:52:45	1:52:50	2-5	G	N
78	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:53:21	1:53:26	2-5	G	T
79	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:54:00	1:54:05	2-5	G	N
80	11	TNT	Neutral	EEG PC	EEG PO	EEG CO	1:54:25	1:54:30	2-5	G	T
81	11	TNT	Neutral	EEG O	EEG P	EEG C	2:14:06	2:14:11	1-pre	97	T
82	11	TNT	Neutral	EEG O	EEG P	EEG C	2:14:25	2:14:30	1-pre	97	N
83	11	TNT	Neutral	EEG O	EEG P	EEG C	2:14:45	2:14:50	1-pre	97	N
84	11	TNT	Neutral	EEG O	EEG P	EEG C	2:15:06	2:15:11	1-pre	97	T
85	11	TNT	Neutral	EEG O	EEG P	EEG C	2:15:28	2:15:33	1-pre	97	T
86	11	TNT	Neutral	EEG O	EEG P	EEG C	2:15:44	2:15:49	1-pre	97	N
87	11	TNT	Neutral	EEG O	EEG P	EEG C	2:16:21	2:16:26	1-pre	97	T
88	11	TNT	Neutral	EEG O	EEG P	EEG C	2:17:10	2:17:15	1-pre	97	N
89	11	TNT	Neutral	EEG O	EEG P	EEG C	2:17:27	2:17:32	1-pre	97	N
90	11	TNT	Neutral	EEG O	EEG P	EEG C	2:17:54	2:17:59	1-pre	97	T

TAPE # 5A

DIGITAL INFORMATION						ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.
		1	2	3	4	5	START	MAKER	Ses. #
91	11	TNT	Neutral	EEG O	EEG P	EEG C	2:18:23	2:18:28	1-pre 97 T
92	11	TNT	Neutral	EEG O	EEG P	EEG C	2:18:45	2:18:50	1-pre 97 N
93	11	TNT	Neutral	EEG O	EEG P	EEG C	2:19:59	2:20:04	1-pre 97 N
94	11	TNT	Neutral	EEG O	EEG P	EEG C	2:20:24	2:20:29	1-pre 97 T
95	11	TNT	Neutral	EEG O	EEG P	EEG C	2:21:02	2:21:07	1-pre 97 N
96	11	TNT	Neutral	EEG O	EEG P	EEG C	2:21:36	2:21:41	1-pre 97 T
97	11	TNT	Neutral	EEG O	EEG P	EEG C	2:21:54	2:21:59	1-pre 97 N
98	11	TNT	Neutral	EEG O	EEG P	EEG C	2:22:20	2:22:25	1-pre 97 T
99	11	TNT	Neutral	EEG O	EEG P	EEG C	2:23:04	2:23:09	1-pre 97 T
100	11	TNT	Neutral	EEG O	EEG P	EEG C	2:23:39	2:23:44	1-pre 97 N
101	11	TNT	Neutral	EEG O	EEG P	EEG C	2:24:39	2:24:44	1-pre 97 T
102	11	TNT	Neutral	EEG O	EEG P	EEG C	2:24:58	2:25:03	1-pre 97 T
103	11	TNT	Neutral	EEG O	EEG P	EEG C	2:25:49	2:25:54	1-pre 97 N
104	11	TNT	Neutral	EEG O	EEG P	EEG C	2:26:11	2:26:16	1-pre 97 N
105	11	TNT	Neutral	EEG O	EEG P	EEG C	2:26:59	2:27:04	1-pre 97 N
106	11	TNT	Neutral	EEG O	EEG P	EEG C	2:27:29	2:27:34	1-pre 97 T
107	11	TNT	Neutral	EEG O	EEG P	EEG C	2:28:19	2:28:24	1-pre 97 N
108	11	TNT	Neutral	EEG O	EEG P	EEG C	2:28:48	2:28:53	1-pre 97 T
109	11	TNT	Neutral	EEG O	EEG P	EEG C	2:29:13	2:29:18	1-pre 97 N
110	11	TNT	Neutral	EEG O	EEG P	EEG C	2:39:01	2:39:06	1-1 97 T
111	11	TNT	Neutral	EEG O	EEG P	EEG C	2:39:23	2:39:28	1-1 97 T
112	11	TNT	Neutral	EEG O	EEG P	EEG C	2:39:45	2:39:50	1-1 97 N
113	11	TNT	Neutral	EEG O	EEG P	EEG C	2:40:15	2:40:20	1-1 97 N
114	11	TNT	Neutral	EEG O	EEG P	EEG C	2:40:51	2:40:56	1-1 97 T
115	11	TNT	Neutral	EEG O	EEG P	EEG C	2:41:37	2:41:42	1-1 97 N
116	11	TNT	Neutral	EEG O	EEG P	EEG C	2:45:36	2:45:51	1-1 97 N
117	11	TNT	Neutral	EEG O	EEG P	EEG C	2:46:10	2:46:15	1-2 97 T
118	11	TNT	Neutral	EEG O	EEG P	EEG C	2:46:45	2:46:50	1-2 97 T
119	11	TNT	Neutral	EEG O	EEG P	EEG C	2:47:18	2:47:23	1-2 97 T
120	11	TNT	Neutral	EEG O	EEG P	EEG C	2:48:11	2:48:16	1-2 97 T

TAPE # 5A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.	
		1	2	3	4	5	START	MAKER	Ses.	#
121	11	TNT	Neutral	EEG O	EEG P	EEG C	2:48:47	2:48:52	1-2	97 N
122	11	TNT	Neutral	EEG O	EEG P	EEG C	2:53:32	2:53:37	1-3	97 N
123	11	TNT	Neutral	EEG O	EEG P	EEG C	2:54:25	2:54:30	1-3	97 T
124	11	TNT	Neutral	EEG O	EEG P	EEG C	2:54:55	2:55:00	1-3	97 T
125	11	TNT	Neutral	EEG O	EEG P	EEG C	2:55:32	2:55:37	1-3	97 N
126	11	TNT	Neutral	EEG O	EEG P	EEG C	2:56:02	2:56:07	1-3	97 T
127	11	TNT	Neutral	EEG O	EEG P	EEG C	2:56:31	2:56:36	1-3	97 N
128	11	TNT	Neutral	EEG O	EEG P	EEG C	3:26:33	3:26:38	3-1	G T
129	11	TNT	Neutral	EEG O	EEG P	EEG C	3:28:16	3:28:21	3-1	G T
130	11	TNT	Neutral	EEG O	EEG P	EEG C	3:29:55	3:30:00	3-1	G N
131	11	TNT	Neutral	EEG O	EEG P	EEG C	3:30:21	3:30:26	3-1	G T
132	11	TNT	Neutral	EEG O	EEG P	EEG C	3:30:53	3:30:58	3-1	G N
133	11	TNT	Neutral	EEG O	EEG P	EEG C	3:35:20	3:35:25	3-2	G N
134	11	TNT	Neutral	EEG O	EEG P	EEG C	3:35:38	3:35:43	3-2	G N
135	11	TNT	Neutral	EEG O	EEG P	EEG C	3:36:10	3:36:15	3-2	G T
136	11	TNT	Neutral	EEG O	EEG P	EEG C	3:36:37	3:36:42	3-2	G T
137	11	TNT	Neutral	EEG O	EEG P	EEG C	3:37:14	3:37:19	3-2	G N
138	11	TNT	Neutral	EEG O	EEG P	EEG C	3:37:45	3:37:50	3-2	G T
139	11	TNT	Neutral	EEG O	EEG P	EEG C	3:39:15	3:39:20	3-3	G N
140	11	TNT	Neutral	EEG O	EEG P	EEG C	3:39:44	3:39:49	3-3	G T
141	11	TNT	Neutral	EEG O	EEG P	EEG C	3:40:06	3:40:11	3-3	G N
142	11	TNT	Neutral	EEG O	EEG P	EEG C	3:40:15	3:40:20	3-3	G T
143	11	TNT	Neutral	EEG O	EEG P	EEG C	3:40:54	3:40:59	3-3	G N
144	11	TNT	Neutral	EEG O	EEG P	EEG C	3:41:15	3:41:20	3-3	G T
145	11	TNT	Neutral	EEG O	EEG P	EEG C	3:42:25	3:42:30	3-3	G N
146	11	TNT	Neutral	EEG O	EEG P	EEG C	3:42:45	3:42:50	3-4	G T
147	11	TNT	Neutral	EEG O	EEG P	EEG C	3:43:15	3:43:20	3-4	G N
148	11	TNT	Neutral	EEG O	EEG P	EEG C	3:43:33	3:43:38	3-4	G T
149	11	TNT	Neutral	EEG O	EEG P	EEG C	3:44:15	3:44:20	3-4	G N
150	11	TNT	Neutral	EEG O	EEG P	EEG C	3:44:35	3:44:40	3-4	G T

TAPE # 5A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.		
		1	2	3	4	5	START	MAKER	Sess	#	
151	11	TNT	Neutral	EEG O	EEG P	EEG C	3:45:47	3:45:52	3-5	G	T
152	11	TNT	Neutral	EEG O	EEG P	EEG C	3:46:07	3:46:12	3-5	G	T
153	11	TNT	Neutral	EEG O	EEG P	EEG C	3:46:48	3:46:53	3-5	G	N
154	11	TNT	Neutral	EEG O	EEG P	EEG C	3:47:15	3:47:20	3-5	G	N
155	11	TNT	Neutral	EEG O	EEG P	EEG C	3:47:30	3:47:35	3-5	G	T
156	11	TNT	Neutral	EEG O	EEG P	EEG C	3:47:54	3:47:59	3-5	G	T
157	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:31:41	4:31:46	4-1	G	T
158	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:32:15	4:32:20	4-1	G	N
159	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:32:46	4:32:51	4-1	G	N
160	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:33:17	4:33:22	4-1	G	T
161	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:33:45	4:33:50	4-1	G	T
162	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:34:09	4:34:14	4-1	G	N
163	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:35:48	4:35:53	4-2	G	N
164	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:36:18	4:36:23	4-2	G	T
165	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:36:50	4:36:55	4-2	G	N
166	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:37:20	4:37:25	4-2	G	T
167	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:37:43	4:37:48	4-2	G	T
168	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:38:06	4:38:11	4-2	G	N
169	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:40:22	4:40:27	4-3	G	T
170	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:40:45	4:40:50	4-3	G	T
171	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:41:15	4:41:20	4-3	G	N
172	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:41:35	4:41:40	4-3	G	N
173	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:41:57	4:42:02	4-3	G	T
174	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:42:29	4:42:34	4-3	G	N
175	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:46:25	4:46:30	4-4	G	N
176	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:47:07	4:47:12	4-4	G	T
177	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:47:49	4:47:54	4-4	G	N
178	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:48:31	4:48:36	4-4	G	T
179	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:49:09	4:49:14	4-4	G	N
180	11	TNT	Neutral	EEG PC	EEG CO	EKG	4:49:48	4:49:53	4-4	G	T

TAPE # 6A

DIGITAL INFORMATION							ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Reccat.	Sess #
		1	2	3	4	5	START	MAKER		
1	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:26:50	-----	BASLINE	
2	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:31:56	00:32:01	2-1	97 T
3	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:32:48	00:32:53	2-1	97 N
4	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:33:20	00:33:25	2-1	97 T
5	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:34:02	00:34:07	2-1	97 N
6	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:34:39	00:34:44	2-1	97 T
7	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:35:18	00:35:23	2-1	97 N
8	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:37:49	00:37:54	2-2	97 T
9	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:38:09	00:38:14	2-2	97 N
10	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:38:41	00:38:46	2-2	97 T
11	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:39:15	00:39:20	2-2	97 N
12	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:39:47	00:39:52	2-2	97 T
13	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:40:32	00:40:37	2-2	97 N
14	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:42:19	00:42:24	2-3	97 T
15	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:42:54	00:42:59	2-3	97 T
16	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:43:26	00:43:31	2-3	97 N
17	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:44:15	00:44:20	2-3	97 T
18	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:44:57	00:45:02	2-3	97 N
19	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	00:45:37	00:45:42	2-3	97 T
20	11	TNT	Neutral	EEG PO	EEG PC	EEG CO	1:06:34	1:06:39	1-pre	100 T
21	11	TNT	Neutral	EEG P	EEG O	EEG C	1:07:05	1:07:12	1-pre	100 N
22	11	TNT	Neutral	EEG P	EEG O	EEG C	1:07:41	1:07:46	1-pre	100 T
23	11	TNT	Neutral	EEG P	EEG O	EEG C	1:09:41	1:09:46	1-pre	100 N
24	11	TNT	Neutral	EEG P	EEG O	EEG C	1:10:27	1:10:32	1-pre	100 T
25	11	TNT	Neutral	EEG P	EEG O	EEG C	1:11:21	1:11:26	1-pre	100 T
26	11	TNT	Neutral	EEG P	EEG O	EEG C	1:11:55	1:12:00	1-pre	100 N
27	11	TNT	Neutral	EEG P	EEG O	EEG C	1:12:15	1:12:20	1-pre	100 T
28	11	TNT	Neutral	EEG P	EEG O	EEG C	1:12:30	1:12:35	1-pre	100 N
29	11	TNT	Neutral	EEG P	EEG O	EEG C	1:12:47	1:12:52	1-pre	100 T
30	11	TNT	Neutral	EEG P	EEG O	EEG C	1:13:03	1:13:08	1-pre	100 N

TAPE # 6A

DIGITAL INFORMATION							ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. Rat.	
		1	2	3	4	5	START	MAKER	Sess #	
31	11	TNT	Neutral	EEG P	EEG O	EEG C	1:13:36	1:13:41	1-pre100	T
32	11	TNT	Neutral	EEG p	EEG O	EEG C	1:13:54	1:13:59	1-pre100	N
33	11	TNT	Neutral	EEG p	EEG O	EEG C	1:14:11	1:14:16	1-pre100	N
34	11	TNT	Neutral	EEG p	EEG O	EEG C	1:14:33	1:14:38	1-pre100	T
35	11	TNT	Neutral	EEG p	EEG O	EEG C	1:15:03	1:15:08	1-pre100	T
36	11	TNT	Neutral	EEG p	EEG O	EEG C	1:15:31	1:15:36	1-pre100	N
37	11	TNT	Neutral	EEG p	EEG O	EEG C	1:16:32	1:16:37	1-pre100	T
38	11	TNT	Neutral	EEG p	EEG O	EEG C	1:16:49	1:16:54	1-pre100	T
39	11	TNT	Neutral	EEG p	EEG O	EEG C	1:17:16	1:17:21	1-pre100	N
40	11	TNT	Neutral	EEG p	EEG O	EEG C	1:17:41	1:17:46	1-pre100	T
41	11	TNT	Neutral	EEG p	EEG O	EEG C	1:18:07	1:18:12	1-pre100	N
42	11	TNT	Neutral	EEG p	EEG O	EEG C	1:18:27	1:18:32	1-pre100	T
43	11	TNT	Neutral	EEG p	EEG O	EEG C	1:19:05	1:19:10	1-pre100	N
44	11	TNT	Neutral	EEG p	EEG O	EEG C	1:19:27	1:19:32	1-pre100	T
45	11	TNT	Neutral	EEG p	EEG O	EEG C	1:19:50	1:19:55	1-pre100	N
46	11	TNT	Neutral	EEG p	EEG O	EEG C	1:20:23	1:20:28	1-pre100	T
47	11	TNT	Neutral	EEG p	EEG O	EEG C	1:20:50	1:20:55	1-pre100	N
48	11	TNT	Neutral	EEG p	EEG O	EEG C	1:22:25	1:22:30	1-1 100	T
49	11	TNT	Neutral	EEG p	EEG O	EEG C	1:23:06	1:23:11	1-1 100	N
50	11	TNT	Neutral	EEG p	EEG O	EEG C	1:23:27	1:23:32	1-1 100	N
51	11	TNT	Neutral	EEG p	EEG O	EEG C	1:23:51	1:23:56	1-1 100	N
52	11	TNT	Neutral	EEG p	EEG O	EEG C	1:24:11	1:24:16	1-1 100	T
53	11	TNT	Neutral	EEG p	EEG O	EEG C	1:24:39	1:24:44	1-1 100	T
54	11	TNT	Neutral	EEG p	EEG O	EEG C	1:26:11	1:26:16	1-2 100	T
55	11	TNT	Neutral	EEG p	EEG O	EEG C	1:26:41	1:26:46	1-2 100	N
56	11	TNT	Neutral	EEG p	EEG O	EEG C	1:27:16	1:27:21	1-2 100	T
57	11	TNT	Neutral	EEG p	EEG O	EEG C	1:27:37	1:27:42	1-2 100	T
58	11	TNT	Neutral	EEG p	EEG O	EEG C	1:28:07	1:28:12	1-2 100	T
59	11	TNT	Neutral	EEG p	EEG O	EEG C	1:28:41	1:28:46	1-2 100	N
60	11	TNT	Neutral	EEG p	EEG O	EEG C	1:30:04	1:30:09	1-3 100	N

TAPE # 6A

DIGITAL INFORMATION						ANALOG INFORMATION				
Line #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec. #	Stat.
		1	2	3	4	5	START	MAKER		
61	11	TNT	Neutral	EEG P	EEG O	EEG C	1:30:49	1:30:54	1-3	100 T
62	11	TNT	Neutral	EEG P	EEG O	EEG C	1:31:46	1:31:51	1-3	100 N
63	11	TNT	Neutral	EEG P	EEG O	EEG C	1:32:16	1:32:21	1-3	100 N
64	11	TNT	Neutral	EEG P	EEG O	EEG C	1:32:45	1:32:50	1-3	100 T
65	11	TNT	Neutral	EEG P	EEG O	EEG C	1:33:16	1:33:21	1-3	100 T
66	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:10:01	2:10:06	2-1	100 T
67	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:10:33	2:10:38	2-1	100 T
68	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:11:06	2:11:11	2-1	100 N
69	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:11:39	2:11:44	2-1	100 N
70	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:12:07	2:12:12	2-1	100 T
71	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:12:34	2:12:39	2-1	100 N
72	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:13:53	2:13:58	2-2	100 T
73	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:14:25	2:14:30	2-2	100 N
74	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:14:58	2:15:03	2-2	100 T
75	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:15:30	2:15:35	2-2	100 N
76	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:16:19	2:16:24	2-2	100 T
77	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:16:57	2:17:02	2-2	100 N
78	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:18:08	2:18:13	2-3	100 T
79	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:18:46	2:18:51	2-3	100 N
80	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:19:28	2:19:33	2-3	100 T
81	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:19:42	2:19:47	2-3	100 T
82	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:20:34	2:20:39	2-3	100 N
83	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:21:02	2:21:07	2-3	100 N
84	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:23:03	2:23:08	2-4	100 T
85	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:23:42	2:23:47	2-4	100 N
86	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:24:13	2:24:18	2-4	100 T
87	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:24:47	2:24:52	2-4	100 T
88	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:25:25	2:25:30	2-4	100 N
89	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:26:06	2:26:11	2-4	100 T
90	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:28:55	2:29:00	2-5	100 T

TAPE # 6A

DIGITAL INFORMATION						ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Rec Rat.	
		1	2	3	4	5	START	MAKER	Ses	#
91	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:29:20	2:29:25	2-5	100 N
92	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:29:45	2:29:50	2-5	100 N
93	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:30:07	2:30:12	2-5	100 T
94	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:30:36	2:30:41	2-5	100 N
95	11	TNT	Neutral	EEG PC	EEG CO	EEG PO	2:31:03	2:31:08	2-5	100 T
96	11	TNT	Neutral	EEG O	EEG C	EEG P	3:08:36	3:08:41	3-1	97 T
97	11	TNT	Neutral	EEG O	EEG C	EEG P	3:09:34	3:09:39	3-1	97 N
98	11	TNT	Neutral	EEG O	EEG C	EEG P	3:10:15	3:10:20	3-1	97 T
99	11	TNT	Neutral	EEG O	EEG C	EEG P	3:10:48	3:10:53	3-1	97 N
100	11	TNT	Neutral	EEG O	EEG C	EEG P	3:11:19	3:11:24	3-1	97 T
101	11	TNT	Neutral	EEG O	EEG C	EEG P	3:11:57	3:12:02	3-1	97 T
102	11	TNT	Neutral	EEG O	EEG C	EEG P	3:14:24	3:14:29	3-2	97 T
103	11	TNT	Neutral	EEG O	EEG C	EEG P	3:14:46	3:14:51	3-2	97 T
104	11	TNT	Neutral	EEG O	EEG C	EEG P	3:15:12	3:15:17	3-2	97 N
105	11	TNT	Neutral	EEG O	EEG C	EEG P	3:15:36	3:15:41	3-2	97 N
106	11	TNT	Neutral	EEG O	EEG C	EEG P	3:16:06	3:16:11	3-2	97 T
107	11	TNT	Neutral	EEG O	EEG C	EEG P	3:16:46	3:16:51	3-2	97 N
108	11	TNT	Neutral	EEG O	EEG C	EEG P	3:19:24	3:19:29	3-3	97 N
109	11	TNT	Neutral	EEG O	EEG C	EEG P	3:19:51	3:19:56	3-3	97 N
110	11	TNT	Neutral	EEG O	EEG C	EEG P	3:20:26	3:20:31	3-3	97 T
111	11	TNT	Neutral	EEG O	EEG C	EEG P	3:20:42	3:20:47	3-3	97 T
112	11	TNT	Neutral	EEG O	EEG C	EEG P	3:21:07	3:21:12	3-3	97 N
113	11	TNT	Neutral	EEG O	EEG C	EEG P	3:21:31	3:21:36	3-3	97 T
114	11	TNT	Neutral	EEG O	EEG C	EEG P	3:53:40	3:53:45	1-pre	103 T
115	11	TNT	Neutral	EEG O	EEG C	EEG P	3:54:10	3:54:15	1-pre	103 N
116	11	TNT	Neutral	EEG O	EEG C	EEG P	3:54:36	3:54:41	1-pre	103 T
117	11	TNT	Neutral	EEG O	EEG C	EEG P	3:54:59	3:55:04	1-pre	103 T
118	11	TNT	Neutral	EEG O	EEG C	EEG P	3:55:33	3:55:38	1-pre	103 N
119	11	TNT	Neutral	EEG O	EEG C	EEG P	3:56:12	3:56:17	1-pre	103 T
120	11	TNT	Neutral	EEG O	EEG C	EEG P	3:56:36	3:56:41	1-pre	103 T

TAPE # 6A

DIGITAL INFORMATION						ANALOG INFORMATION			
File #	# of Rec.	CONTENTS IN-CHANNEL					TIME CODE		Rec. Rat.
		1	2	3	4	5	START	MAKER	Ses #
121	11	TNT	Neutral	EEG O	EEG C	EEG P	3:57:03	3:57:08	1-pre103 N
122	11	TNT	Neutral	EEG O	EEG C	EEG P	3:57:38	3:57:43	1-pre103 T
123	11	TNT	Neutral	EEG O	EEG C	EEG P	3:58:33	3:58:38	1-pre103 T
124	11	TNT	Neutral	EEG O	EEG C	EEG P	3:59:03	3:59:08	1-pre103 N
125	11	TNT	Neutral	EEG O	EEG C	EEG P	3:59:41	3:59:46	1-pre103 T
126	11	TNT	Neutral	EEG O	EEG C	EEG P	3:59:59	4:00:04	1-pre103 T
127	11	TNT	Neutral	EEG O	EEG C	EEG P	4:00:13	4:00:18	1-pre103 N
128	11	TNT	Neutral	EEG O	EEG C	EEG P	4:00:40	4:00:45	1-pre103 T
129	11	TNT	Neutral	EEG O	EEG C	EEG P	4:01:20	4:01:25	1-pre103 N
130	11	TNT	Neutral	EEG O	EEG C	EEG P	4:01:54	4:01:59	1-pre103 T
131	11	TNT	Neutral	EEG O	EEG C	EEG P	4:02:16	4:02:21	1-pre103 N
132	11	TNT	Neutral	EEG O	EEG C	EEG P	4:02:42	4:02:47	1-pre103 N
133	11	TNT	Neutral	EEG O	EEG C	EEG P	4:04:01	4:04:06	1-pre103 T
134	11	TNT	Neutral	EEG O	EEG C	EEG P	4:04:22	4:04:27	1-pre103 N
135	11	TNT	Neutral	EEG O	EEG C	EEG P	4:04:36	4:04:41	1-pre103 N
136	11	TNT	Neutral	EEG O	EEG C	EEG P	4:05:02	4:05:07	1-pre103 T
137	11	TNT	Neutral	EEG O	EEG C	EEG P	4:05:22	4:05:27	1-pre103 T
138	11	TNT	Neutral	EEG O	EEG C	EEG P	4:05:41	4:05:46	1-pre103 T
139	11	TNT	Neutral	EEG O	EEG C	EEG P	4:06:01	4:06:06	1-pre103 N
140	11	TNT	Neutral	EEG O	EEG C	EEG P	4:06:34	4:06:39	1-pre103 T
141	11	TNT	Neutral	EEG O	EEG C	EEG P	4:07:13	4:07:18	1-pre103 N
142	11	TNT	Neutral	EEG O	EEG C	EEG P	4:07:40	4:07:45	1-pre103 T
143	11	TNT	Neutral	EEG O	EEG C	EEG P	4:10:31	4:10:36	1-1 103 T
144	11	TNT	Neutral	EEG O	EEG C	EEG P	4:10:51	4:10:56	1-1 103 T
145	11	TNT	Neutral	EEG O	EEG C	EEG P	4:11:25	4:11:30	1-1 103 N
146	11	TNT	Neutral	EEG O	EEG C	EEG P	4:11:41	4:11:46	1-1 103 N
147	11	TNT	Neutral	EEG O	EEG C	EEG P	4:12:03	4:12:08	1-1 103 T
148	11	TNT	Neutral	EEG O	EEG C	EEG P	4:12:26	4:12:31	1-1 103 N
149	11	TNT	Neutral	EEG O	EEG C	EEG P	4:13:59	4:14:04	1-2 103 T
150	11	TNT	Neutral	EEG O	EEG C	EEG P	4:14:29	4:14:34	1-2 103 N

TAPE # 6A

DIGITAL INFORMATION							ANALOG INFORMATION				
File #	# of Rec.	CONTENTS IN CHANNEL					TIME CODE		Recorder	Tape #	Status
		1	2	3	4	5	START	MAKER			
151	11	TNT	Neutral	EEG O	EEG C	EEG P	4:14:55	4:15:00	1-2	103	T
152	11	TNT	Neutral	EEG O	EEG C	EEG P	4:15:17	4:15:22	1-2	103	N
153	11	TNT	Neutral	EEG O	EEG C	EEG P	4:15:43	4:15:48	1-2	103	T
154	11	TNT	Neutral	EEG O	EEG C	EEG P	4:16:12	4:16:17	1-2	103	N
155	11	TNT	Neutral	EEG O	EEG C	EEG P	4:17:47	4:17:52	1-5	103	N
156	11	TNT	Neutral	EEG O	EEG C	EEG P	4:18:15	4:18:20	1-3	103	N
157	11	TNT	Neutral	EEG O	EEG C	EEG P	4:19:02	4:19:07	1-3	103	T
158	11	TNT	Neutral	EEG O	EEG C	EEG P	4:19:41	4:19:46	1-3	103	N
159	11	TNT	Neutral	EEG O	EEG C	EEG P	4:20:00	4:20:05	1-5	103	T
160	11	TNT	Neutral	EEG O	EEG C	EEG P	4:20:34	4:20:39	1-3	103	T
161	11	TNT	Neutral	EEG O	EEG C	EEG P	4:23:15	4:23:20	1-4	103	N
162	11	TNT	Neutral	EEG O	EEG C	EEG P	4:23:46	4:23:51	1-4	103	N
163	11	TNT	Neutral	EEG O	EEG C	EEG P	4:24:39	4:24:44	1-4	103	T
164	11	TNT	Neutral	EEG O	EEG C	EEG P	4:25:19	4:25:24	1-4	103	T
165	11	TNT	Neutral	EEG O	EEG C	EEG P	4:25:57	4:26:02	1-4	103	N
166	11	TNT	Neutral	EEG O	EEG C	EEG P	4:26:35	4:26:40	1-4	103	T
167	11	TNT	Neutral	EEG O	EEG C	EEG P	4:27:48	4:27:53	1-5	103	N
168	11	TNT	Neutral	EEG O	EEG C	EEG P	4:28:15	4:28:20	1-5	103	N
169	11	TNT	Neutral	EEG O	EEG C	EEG P	4:28:42	4:28:47	1-5	103	T
170	11	TNT	Neutral	EEG O	EEG C	EEG P	4:29:03	4:29:08	1-5	103	T
171	11	TNT	Neutral	EEG O	EEG C	EEG P	4:29:18	4:29:23	1-5	103	N
172	11	TNT	Neutral	EEG O	EEG C	EEG P	4:29:55	4:30:00	1-5	103	T


TASK IV objectives were to:

1. Supply sufficient documentation detailing procedures needed for data extraction from the digital tapes. This information pertains to tape formatting and unpacking procedures, that are required, and have to be followed for magnetic tape dumps, processing and so on.


The following steps are needed to unpack and de-multiplex the 12-bit data into a 16-bit word data for ease of handling.

- 1) Get runtime parameters:
 - a) Number of channels recorded = NCH
 - b) Number of data valves per channel = NWDCH
- 2) Calculate number of bytes per input record according to runtime parameters.
 - a) $NBUT = (NCH * NWDCH) * 3/2$
- 3) Unpack and de-multiplex each record in accordance to these formats:

Unpack:



De-multiplex:



0001 PROGRAM UNPACK.

GENERATE TAPE FILE OF DATA COLLECTED BY THE DIGI-DATA
RERECORDER FOR SANDIA AND SENSOR ARRAY PROJECTS.

```

0002 LOGICAL*1 BYTEA(2), ISD
0003 INTCR OFILE,JARY(4076), JARY(4076), ISTAT(2)
0004 INTEGER IN(256,B),OUT(256,B)
0005 EQUIVALENCE (JARY(1),BYTEA(1)),(ISTAT(1),ISD)
0006 EQUIVALENCE (JARY(1),IN(1,B))
0007 DATA LONG/117,LUH/107,MAXLEN/192/

```

GET FINAL PARAMETERS

```
0008 ACCEPT 100,NCH,INCHL,MODE,NRTD,IFILE,NSKP,LFILE,OFILF
0009 FORMAT(815)
0010 IFILE=IFILE
```

0011 PRINT IG1,NCH,NDCN1,PCDE,IFILL,NSKP

0012 101 FCBM (7-59, 'DIGI-DATA WRAPPER', //

STIQ, '4 CF CHANNELS RECORDED', T50, 15, /,

DETERMINE # OF BYTES (FRAMES) PER INPUT RECORD

6015 $\text{Ni}=\text{C}_2\text{H}_4\text{V}_2\text{O}_2\text{CH}$

0016 0047:040

6517
 100-111042

8100 IF(MODE, EQ, 2)NB=NH*2/2

0019
PRINT 102, MS. MICH. NDS

0020 102 FORMAT(T10, '# OF CYLES (FRAMES) PER INPUT RECORD', T50.15. /,

C DATA RECORDS TOO LARGE

0022 103. M. 4. 110

```

00023      FORMAT(' EXCEEDED ARRAY SIZE OF', I5, ' BYTES')

```

0024 5103

INITIAL TALE POSITION

0025 S CALL INTERPOL UNIT, 194, 1ST ST)

6026 CALL MTRNGS (LUND, 1981, 151AT)

C . . SKIP HEADER RECORD FOR THIS FILE

0027 CALL MTCOFF(CUR), FILE=-1, ISW, ISJAT)

0029 IF(IEB LT 0)GC TO 99

UNPACK.FTN

/11. DLOCKS/WR

CALL MTKPPF(LUND, OFILE-1, ISW, ISTAT)

IF (ISB.LT.0) GO TO 70

HOW MANY INPUT RECORDS TO SKIP

CALL MTKPPR(LUNI, NCKP, ISW, ISTAT)

GET DATA IN(CORD

IR=0

CONTINUE

JX=1

CALL MTREAD(LUNI, DYIEA, NB, ISW, ISTAT)

IF (ISD.EQ.-10) GO TO 30

IF (ISB.EQ.-11) GO TO 30

IF (ISD.LT.0) GO TO 91

UNPACK DIGIT DATA FORMAT

CALL UNPKDD(DYIEA, NB, IARY, MODE)

SE-MULTIPLY X

CALL DMUXDD(IARY, JARY, NCH, NUDDCH)

OUTPUT TO TAPE

DO 33 IC=1, NDATA, NUDDCH

CALL MTWRIT(LUND, JARY(IC), NDB, ISW, ISTAT)

IF (ISD.LT.0) GO TO 92

CONTINUE

CHECK FOR # RECORDS/FILE

IR=IR+1

IF (IR.LT.NRTD) GO TO 12

CALL MTKPPF(LUNI, 1, ISW, ISTAT)

END OF DATA THIS FILE

PRINT 106, IR

FORMAT(10, 'TOTAL RECORDS TO TAPE', T50, 15)

CALL MTKGOF(LUND, ISW, ISTAT)

IF IL=IFIL+1

IN = 0

IF (IFIL.LE.LFILE) GO TO 12

CALL MTKGOF(LUND, ISW, ISTAT)

CALL MTKGOF(LUNI, ISW, ISTAT)

CALL MTKGOF(LUND, ISW, ISTAT)

STOP

ERROR HANDLERS

PRINT 700, ISE

FORMAT(' TAPE ERROR ON HEADER READ', 15)

STOP

FORTRAN IV-PLUS V02-51

UNPACK.FTN /TR:BLOCKS/WR

12:32:30

24-AUG-79

PAGE 3

```
0061 91 PRINT '10,15B,IR
0062 910 FORMAT(' TAPE ERROR ON DATA READ',I5,10X,'RECORD',I5)
0063 GO TO 18
0064 92 PRINT '20,15B,IR
0065 920 FORMAT(' TAPE ERROR ON DATA WRITE',I5,10X,'RECORD',I5)
0066 C
0067 END
```


PROGRAM SECTIONS

NUMBER	NAME	SIZE	ATTRIBUTES
1	\$CODE1	001416 391	RW, I, CON, LCL
2	\$IDATA	000004 2	RW, D, CON, LCL
3	\$IDATA	001272 349	RW, D, CON, LCL
4	\$VARS	030052 10261	RW, D, CON, LCL
5	\$TEMPS	000002 1	RW, D, CON, LCL

VARIABLES

NAME	TYPE	ADDRESS	NAME	TYPE	ADDRESS	NAME	TYPE	ADDRESS	NAME	TYPE	ADDRESS
IC	I*2	4-050050	IFIL	I*2	4-050032	IR	I*2	4-050024	ISB	L*1	4-0253
ISW	I*2	4-050042	JX	I*2	4-050046	LUNI	I*2	4-050030	LUND	I*2	4-0503
MAXND	I*2	4-050012	MODE	I*2	4-050020	NCH	I*2	4-050034	NDATA	I*2	4-0503
NDJ	I*2	4-050040	NRTD	I*2	4-050022	NWDCH	I*2	4-050026	OFFILE	I*2	4-0203

ARRAYS

NAME	TYPE	ADDRESS	SIZE	DIMENSIONS
BYTEA	L*1	4-000000	000002	1 (2)
IARY	I*2	4-070006	020000	4096 (4096)
IN	I*2	4-000000	010000	2048 (256,8)
ISTAT	I*2	4-030000	000004	2 (2)
JARY	I*2	4-000000	020000	4096 (4096)
OUT	I*2	4-040006	010000	2048 (256,8)

LABELS

LABEL	ADDRESS	LABEL	ADDRESS	LABEL	ADDRESS	LABEL	ADDRESS
5	1-000512	12	1-000654	18	1-000726	19	1-001062
33	1-001236	90	1-001236	71	1-001300	92	1-001350
101	3-000004	102	3-000440	103	3-000636	105	3-000276
900	3-000742	910	3-001002	920	3-001054		

FUNCTIONS AND SUBROUTINES REFERENCED

DMUXDD MTRWDD MTRWDD MTSKPF MTSKPF MTRWDD MTRWDD UNPACKD

TOTAL SPACE ALLOCATED = 052770 11004

NO FPP INSTRUCTIONS GENERATED

SYO: [234, 222]UNPACK, SYO: [234, 222]UNPACK/-SP=SYO: [234, 222]UNPACK